

SofCare® HeelCare® Boot Series



Convenient. Complete. Effective.

Stryker's SofCare family of products allows for a convenient, complete and effective way to enhance pressure redistribution and patient comfort while improving hospital outcomes and helping to prevent facility acquired pressure ulcers.

Stryker has adapted its SofCare air-channeling technology into light weight, comfortable and easy to use therapeutic air boots for the heel and foot. The adjustable SofCare HeelCare Boot series lifts the heel off the bed surface, helping to prevent pressure ulcers while allowing air circulation.

Open Design

Provides added air circulation for a cooler, comfortable environment

Adjustable Straps

Quick and easy set up while customizing fit for each patient

Nylon Exterior

Easily cleaned and helps allow resistance free movement reducing sheering effects

Soft Scalloped Edges

Provide added comfort and no sharp edges against delicate skin

Zero Heel Pressure

Elevates heel totally off the surface for complete pressure relief

Preinflated

No need to reinflate due to urethane coated materials that fully seal the heel boot and minimize stretch



Elite Heel Boot

HCB3000

- Scalloped heel design that protects against pressure with movement
- Upgraded material
- Adjustable straps for customized fit



Universal Heel Boot

HCB2001

- Designed to fit most foot sizes
- Easy access and heel visibility
- Adjustable straps for customized fit




Long Heel Boot

HCB203SL

- Longer length for more complete lower extremity protection
- Traction grommet for orthopedic application
- Fits calf 12 to 14 inches

Air-Channeling Technology

SofCare air-channelling technology redistributes tissue interface pressure to aid in prevention and treatment of pressure ulcers. The channel design allows for air circulation helping to manage patient microclimate.



SofCare Air-Channeling Technology

SofCare HeelCare Boot Series



SofCare HeelCare Boot Series

- Channeled tubular support
- Soft inner material and scalloped edges
- Easy access and heel visibility
- Adjustable straps for customized fit
- Zero heel Pressure

Specifications

Boot	Elite HeelCare Boot	Universal Heel	Long Heel Boot
Model Number	HCB3000	HCB2001	HCB204ML
Weight	4 oz	4 oz	5 oz
Packaging	6/carton	6/carton	6/carton
Materials	Nylon	Poly Nonwoven	Poly Nonwoven
Warranty	30 days	30 days	30 days

Single Patient Use Only

Ask about our other SofCare 3-layer air-channeling technology products. Stryker reserves the right to change specifications without notice.

3800 E. Centre Ave.
Portage, MI 49002 U.S.A.
t: 269 389 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com