

# Spirit<sup>®</sup>

## Behavioral Health Platform Bed



### Contemporary design with robust construction

Our behavioral health platform bed combines contemporary design with robust construction. Tamper-resistant hardware, durable frame materials and a large weight capacity make this bed ideal for behavioral health facilities requiring a resilient bed frame. The flat surfaces of the bed make it easy to clean and sanitize.

#### Headboard and footboard are high density polyethylene (HDPE)

- Easy to clean; stain and odor resistant
- Will not chip, swell, rot or crack

#### Tamper-resistant fasteners used throughout bed

#### Standard model (80")

- Restraint loops (optional)
- Rubber feet
- 3-inch headboard / 3-inch footboard
- Bottom pan

#### Extended model (84")

- Restraint loops
- Rubber feet
- 7-inch headboard / 3-inch footboard
- Bottom pan

Specifications	Standard model	Extended model
<b>Model number</b>	5800	5800
<b>Overall length</b>	81.5" (207 cm)	85.5" (217.2 cm)
<b>Mattress length</b>	80" (203.2 cm)	84" (218.3 cm)
<b>Overall width</b>	36" (91.4 cm)	36" (91.4 cm)
<b>Weight capacity</b>	1000 lb. (454 kg)	1000 lb. (454 kg)
<b>Height (to liter top)</b>	16" (40.6 cm)	16" (40.6 cm)

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Spirit, Stryker. All other trademarks are trademarks of their respective owners or holder.

3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770

[stryker.com](http://stryker.com)

Mkt Lit-1049 08 JAN 2015 Rev E  
Copyright © 2018 Stryker  
Printed in U.S.A