

# Spirit<sup>®</sup> Behavioral Health Bed



## Low height, high security.

Behavioral health departments can place unique demands on hospital employees. We worked to create a solution with features that help to satisfy these demands and increase security within these specific units when partnered with specialized hospital processes. In addition to the specialized security features that can be added to Spirit beds, all beds have a 10.75-inch low bed height. This low height helps to reduce fall injuries and increase patient comfort to create a more secure sleeping environment.



Any or all patient controls can be locked to restrict use with the patient lockout buttons. Behavioral health beds are equipped with an additional security code lock-out.



Optional bed exit alarm monitors patient's bed exit activity in single-zone or three zone.



Nonremoveable mattress deck.

**Siderails** comply with the FDA's monitored Entrapment Zones 1 – 4

**Filled siderails, headboard and footboard** to help prevent harm

**Patient lockout** restricts all bed functions

**Nonremoveable** mattress deck, headboard, and footboard

**Tamper-resistant hardware**

**Permanent control box cover**

**Removable brake pedal**

**Shortened power cord**

# Spirit Behavioral Health Bed



## Spirit Plus Behavioral Health Bed

- Battery back-up
- Patient lock-outs
- Dual lock-and-steer brake system
- Footboard control panel
- Auto contour
- Low chair positioning
- Trendelenburg/reverse Trendelenburg
- Total bed lock-out
- Molded plastic headboard/footboard
- Embedded patient controls in siderails
- Integrated nurse call button
- Power status indicator
- High quad rails
- Filled rails and boards
- 4-inch tool-less bed extension
- Behavioral health package (tamper-resistant hardware, nonremovable headboard and footboard, removable brake pedal)
- Security lockout
- Shortened power cord

## Spirit Select Base Behavioral Health Bed

includes all features listed above, plus:

- Backlit staff controls
- Level-all bed reset
- Built-in scale
- Integrated single-zone bed exit alarm (optional 3-zone alarm)
- Programmable night light

## Spirit Select Loaded Behavioral Health Bed

includes all features listed above, plus:

- Watchdog Bed Status Indicator
- Obstruction sensors
- Nurse reminder
- Additional colored night lights
- 120V AC accessory outlet

## Specifications

<b>Model number</b>	5600/5700 (plus/select)
<b>Overall length</b>	91.75" (233 cm)
<b>Overall width</b>	40" (101.6 cm)
<b>Safe working load</b>	500 lb. (227 kg)
<b>Height range (to litter top)</b>	
High	34" (86 cm)
Low	10.75" (27 cm)
<b>Litter positioning</b>	
Knee gatch	0°–25°
Trend./Reverse Trend.	14° (+/-)
Backrest	0°–60°
Foot section	0°–10°
<b>Patient surface</b>	
Minimum	35" x 80" (89 x 203 cm)
Maximum	35" x 84" (89 x 213 cm)
<b>Caster diameter</b>	5" (12.7 cm)

## Electronics

<b>Volts</b>	120 VAC
--------------	---------

3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770

[stryker.com](http://stryker.com)

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Spirit, Stryker, Watchdog. All other trademarks are trademarks of their respective owners or holder.