

# Spirit Select<sup>®</sup> Postpartum Bed



## Safety focused for mom and baby

In conjunction with effective protocols, the Spirit Select Postpartum Bed can enhance the care and safety of mothers and babies throughout their time in the hospital.

The postpartum siderail allows patients to control the height of the bed, as well as the position of the siderails from the inside control panel located on the interior of the upper siderails.

Open architecture, easy-to-clean bed frame, and infinitely adjustable positions can help caregivers tend to their patients' needs while reducing risk to their own safety.



Having the ability to lower the siderails easily with the touch of a button helps give new mothers the independence and security they desire when transferring their babies between the bed and the bassinet before and after feedings.



Bed exit alarm helps protect staff and patients from common alarm hazards such as false alarms, alarm desensitization and nonrestoration of settings.



Any or all patient controls can be locked to restrict use with the patient lockout buttons.

**10.75-inch low bed height** allows patient to enter and exit with stability

**Programmable night light** adds an additional element of safety during patient transfers

**Built-in smart scale** allows for accuracy\* and is equipped with a weight-change alarm

**Nurse call connectivity** provides better staff and patient communication

\*  $\pm 3\%$  for patients who weigh 110 to 500 lb.  
 $\pm 3$  lb. for patients who weigh less than 110 lb.

# Spirit Select Postpartum Bed



## Standard features

- Battery back-up
- Patient lock-outs
- Dual lock-and-steer brake system
- Footboard control panel
- Auto contour
- Low chair positioning
- Trendelenburg/reverse Trendelenburg
- Emergency CPR — automatic
- Emergency CPR — manual
- Total bed lock-out
- Molded plastic headboard/footboard
- Embedded patient controls in siderails
- Integrated nurse call button
- Backlit staff controls
- Power status indicator
- Level-all bed reset
- Built-in scale
- Integrated single-zone bed exit alarm
- Programmable night light
- 4-inch tool-less bed extension (option)
- 120 V AC accessory outlet
- Obstruction sensors (option)
- High quad rails with filled foot rail
- Patient controlled foot rail release
- Patient controlled bed height

## Specifications

<b>Model number</b>	5700
<b>Overall length</b>	91.75" (233 cm)
<b>Overall width</b>	40" (101.6 cm)
<b>Weight capacity</b>	500 lb. (227 kg)
<b>Height range (to litter top)</b>	
High	34" (86.4 cm)
Low	10.75" (27.3 cm)
<b>Litter positioning</b>	
Knee gatch	0°-25°
Trend./Reverse Trend.	14° (+/-)
Backrest	0°-60°
Foot section	0°-10°
<b>Patient surface</b>	
Minimum	35" x 80" (89 x 203 cm)
Maximum	35" x 84" (89 x 213 cm)
<b>Caster diameter</b>	5" (12.7 cm)

## Electronics

<b>Volts</b>	120 VAC
--------------	---------

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Spirit Select, Stryker, Watchdog. All other trademarks are trademarks of their respective owners or holder.

3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770

[patientcare.stryker.com](http://patientcare.stryker.com)

Mkt Lit-1051 08 JAN 2015 Rev C.3  
Copyright © 2017 Stryker  
Printed in U.S.A