

Spirit Select[®] Specialty Bed



High confidence, low height

With the right programs and processes in place, the Spirit Select Bed is designed to help keep both patients and staff safe. Spirit Select is a specialty low acute care bed with a low bed height of 10.75 inches from floor to litter deck. Low height enhances ingress and egress from the bed by adding additional stability. The Spirit Select bed can also help improve staff efficiency and patient outcomes with nurse call connectivity, fall alarm, 32-character LED display, battery backup and Watchdog Bed Status Indicator.



Short hand rail allows patient support during transfers.



Bed exit alarm helps protect staff and patients from common alarm hazards such as false alarms, alarm desensitization and nonrestoration of settings.



Any or all patient controls can be locked to restrict use with the patient lockout buttons.

10.75-inch low bed height allows patient to enter and exit with stability

Programmable night light adds an additional element of safety during patient transfers

Built-in smart scale allows for accuracy* and is equipped with a weight-change alarm

Nurse call connectivity provides better staff and patient communication

* $\pm 3\%$ for patients who weigh 110 to 500 lb.
 ± 3 lb. for patients who weigh less than 110 lb.

Spirit Select Specialty Bed



Standard features

- Battery back-up
- Patient lock-outs
- Dual lock-and-steer brake system
- Footboard control panel
- Auto contour
- Low chair positioning
- Trendelenburg/reverse Trendelenburg
- Emergency CPR — automatic
- Emergency CPR — manual
- High quad rails
- Total bed lock-out
- Molded plastic headboard/footboard
- Embedded patient controls in siderails
- Integrated nurse call button
- Backlit staff controls
- Power status indicator
- Level-all bed reset
- Built-in scale
- Integrated single-zone bed exit alarm
- Programmable night light

Optional features

- 120 VAC accessory outlet
- Obstruction sensors
- Watchdog Bed Status Indicator
- 4-inch tool-less bed extension

Specifications

Model number	5700
Overall length	91.75" (233 cm)
Overall width	40" (101.6 cm)
Safe working load	500 lb. (227 kg)
Height range (to litter top)	
High	34" (86 cm)
Low	10.75" (27 cm)
Litter positioning	
Knee gatch	0°-25°
Trend./Reverse Trend.	14° (+/-)
Backrest	0°-60°
Foot section	0°-10°
Patient surface	
Minimum	35" x 80" (89 x 203 cm)
Maximum	35" x 84" (89 x 213 cm)
Caster diameter	5" (12.7 cm)

Electronics

Volts	120 VAC
--------------	---------

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Spirit Select, Stryker, Watchdog. All other trademarks are trademarks of their respective owners or holder.

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

patientcare.stryker.com

Mkt Lit-1142 09 JUN 2015 Rev B.3
Copyright © 2017 Stryker
Printed in U.S.A