

# Symmetry® II

## Patient and Guest Seating



Highback Chair

Lowback Chair With Arms

Lowback Chair Without Arms

**High-density Ultracell foam provides lasting shape and support**

**Fully-extended arms permit easy entry and exit**

**Accentuated lumbar design supports lower back**

**Open design facilitates cleaning**

**Wall saver design protects integrity of chair and walls**

**Reinforced joints provide increased stability and strength**

**Limited lifetime warranty**

### Comprehensive Patient Room Seating Solutions

Symmetry seating offers comfortable, stylish alternatives for your patient rooms. Each chair's ergonomic design provides lumbar, leg and hip support while the highly-resilient construction complements clinical standards.



Flexback Highback Chair



Steel Chair



Ottoman

# Symmetry® II

## Patient and Guest Seating

**stryker®**



Highback Chair

Flexback Highback Chair

Lowback Chair  
With Arms

Lowback Chair  
Without Arms

Steel Chair

### Medical

3800 E. Centre Ave.  
Portage, MI 49002  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

www.stryker.com

### Standard Features

- Extended arm caps make chair easy to grasp and facilitates patient exit and entry
- Accentuated lumbar design supports lower back
- Flexible back design allows for increased comfort (Flexback Highback chair only)
- Meets flammability standards of CAL 117
- Meets ANSI/BIFMA X5.1-2002 per tests conducted in a certified test facility

### Standard Features (Steel Modular Chair)

- Strong, lightweight tube steel frame
- Durable plastic arm rests
- Stackable design
- Meets flammability standards of CAL 117
- Meets ANSI/BIFMA X5.1-2002 per tests conducted in a certified test facility

### Specifications

• Style	<b>Highback</b> Part #3500-000-110 (18") Part #3500-000-130 (21")	<b>Flexback Highback</b> Part #3500-000-320 (18") Part #3500-000-340 (21")
• Overall Width	24" (61 cm)	24" (61 cm)
• Overall Height	41"/44" (103/110 cm)	40.5"/43.5" (103/111 cm)
• Depth		
Overall	26" (66 cm)	26" (66 cm)
Seat	17.5" (45 cm)	18" (46 cm)
• Weight	21/24 lbs (10/11 kg)	22/24 lbs (10/11 kg)
• Seat Width	21" (53 cm)	21" (53 cm)
• Back Height From Seat	23" (58 cm)	23" (58 cm)
• Seat Height From Floor	18"/21" (45/53 cm)	18" (46 cm)
• Arm Height From Floor	25.5"/28.5" (65/74 cm)	25.5"/28.5" (65/72 cm)

• Style	<b>Lowback w/Arms</b> Part #3500-000-150 (18") Part #3500-000-170 (21")	<b>Lowback w/o Arms</b> Part #3500-000-220	<b>Steel Chair</b> Part #3500-000-520
• Overall Width	24" (61 cm)	22.5" (57 cm)	23" (58 cm)
• Overall Height	35"/38" (88/96 cm)	35" (88 cm)	32" (81 cm)
• Depth			
Overall	26" (66 cm)	26" (66 cm)	20" (51 cm)
Seat	18" (46 cm)	18" (46 cm)	17.5" (45 cm)
• Weight	21/22 lbs (10 kg)	19 lbs (9 kg)	14 lbs (6 kg)
• Seat Width	21" (53 cm)	21" (53 cm)	18" (46 cm)
• Back Height From Seat	17" (43 cm)	17" (43 cm)	14" (36 cm)
• Seat Height From Floor	18"/21" (45/53 cm)	18" (45 cm)	18" (46 cm)
• Arm Height From Floor	25.5"/28.5" (65/74 cm)	N/A	26" (66 cm)

Stryker reserves the right to change specifications without notice.  
Furniture is manufactured exclusively for Stryker Medical.

### Warranty

- Limited lifetime warranty
- Warrantied 300 lbs (136 kg) weight capacity

### Optional Features

- Choose from a variety of stains and fabrics
- CAL 133