

Symmetry® Plus Patient Room Recliners

Versatile recliners for a healing environment.



Stryker's stylish Patient Room Recliner **maximizes patient and guest comfort**, while providing uncompromised durability. Choose from Rocker or Wall-Saver models, with open or upholstered arms available on each model.



Step-Lock Brakes are easy to access to help improve caregiver efficiency.



Rear Push Bar makes transporting a patient easier.



Rolling footrest provides added comfort while in transit.

Backrest stops in infinite number of recline positions

Posture-correct lumbar and lateral cushioning

Pull-out foot tray gives patients added comfort while in transit

Rocker or Wall-Saver models available

Extended, wide arm caps facilitate easy entry and exit

Symmetry Plus Patient Room Recliners



Standard Features

- Stationary or mobile
- Mobile with 3 inch front rigid/rear swivel locking casters
- Mid-ottoman cushion provides extra support when feet are elevated
- Chair construction includes metal seat frame and sinuous wire springs
- Highly resilient polyurethane foam cushions
- Removable back cushion enables thorough cleaning and service
- Pull-out foot tray gives patients added comfort while in transit
- Rigid push bar for easy maneuverability
- Stationary arms and fixed aides provide stability for drainage catheters or electromechanical mounts (Open-Arm only)
- Egress assist polyurethane arm caps
- Meets flammability standards of BIFMA F1-1978 — revised 1980 and CAL 117
- Meets ANSI//BIFMA X5.4 2005 per tests conducted in a certified test facility

Optional Features

- Open arm or upholstered arm
- Egress assist wood arm caps
- Hinged push bar for maneuverability
- Four-wheel swivel with locking casters
- Choose from a variety of stains and fabrics
- CAL 133 in select fabrics
- Foley Bag hooks (Wall-Saver only)
- 3 inch front/rear swivel locking casters

Specifications

| | Rocker Recliner | Wall-Saver |
|------------------------|------------------------|----------------------|
| Part Number | 3500-000-750 | 3500-000-730 |
| Overall Width | 31" (78.7 cm) | 31" (78.7 cm) |
| Overall Height | 45.5" (115.6 cm) | 47.5" (120.7 cm) |
| Depth | | |
| Overall | 41.5" (105.4 cm) | 40.5" (102.9 cm) |
| Seat | 19.5" (49.5 cm) | 19.5" (49.5 cm) |
| Recline / Leg Rest Up | 70" (177.8 cm) | 70" (177.8 cm) |
| Weight | 138.2 lb. (62.7 kg) | 130.5 lb. (331.5 kg) |
| Seat Width | 23" (58.4 cm) | 23" (58.4 cm) |
| Back Height From Seat | 32.5" (82.6 cm) | 33.5" (85.1 cm) |
| Seat Height From Floor | 20.5" (52.1 cm) | 20.5" (52.1 cm) |
| Arm Height From Floor | 27" (68.6 cm) | 27" (68.6 cm) |

Warranty

Limited lifetime warranty
Warrantied 350 lb. (158.8 kg) weight capacity

Rocker Recliner

Rear lock lever allows chair to remain stationary

Wall-Saver

Reclines away from the wall to save space

*Stryker reserves the right to change specifications without notice.
Furniture is manufactured exclusively for Stryker Medical.*