

# Tru-Fit® Overbed Tables



Split-top with storage compartment  
(available on split-top or single-top)

The sleek design of the award-winning Tru-Fit Overbed Table easily slides under beds, stretchers and mobile recliners. Available in several laminate and base colors, The Stryker single-top and split-top overbed tables coordinate with any room decor.



Single-top with vanity



Single-top without vanity



Split-top overbed table

**Low-profile, U-shaped base** provides maximum versatility.

**Infinite adjustments** between low and high heights

**High-impact polystyrene cover** protects the base and conceals the caster

**Dura-edge™** protects the table and contains spills

**Available with** vanity compartment and acrylic mirror

**Storage compartment option** provides additional storage space

# Tru-Fit

## Overbed Tables

stryker



Single-top without vanity



Single-top with vanity



Split-top overbed table with vanity



Split-top with storage compartment

### Standard features

- 1.46 inch (3.7 cm) diameter twin swivel casters
- Durable high-impact polystyrene base cover
- Height adjustment to 43.75 inches (111 cm) (single-top only)
- Height adjustment to 45.5 inches (115 cm) (split-top only)
- Infinite adjustments between low and high heights
- Anodized aluminium column
- Tabletop with Dura-Edge
- Secondary laminate top (split-top only)
- Storage vanity with acrylic mirror

### Optional features

- Choose from six standard laminates and four standard base/Dura-Edge colors
- Storage compartment

### Specifications

	Single-top without vanity	Single-top with vanity	Split-top with vanity
<b>Part #</b>	3150-000-100	3150-000-200	3150-000-300
<b>Height</b>			
<b>High</b>	43.75" (111.1 cm)	43.75" (111.1 cm)	45.5" (115.6 cm)
<b>Low</b>	27" (68.6 cm)	27" (68.6 cm)	29" (73.7 cm)
<b>Weight capacity</b>	75 lb. (34 kg)	75 lb. (34 kg)	75 lb. (34 kg)
<b>Weight</b>	43 lb. (20 kg)	50 lb. (23 kg)	61 lb. (28 kg)
<b>Overall table top</b>			
<b>Length</b>	31" (79 cm)	31" (79 cm)	31" (79 cm)
<b>Width</b>	17.5" (44.5 cm)	17.5" (44.5 cm)	17.5" (44.5 cm)
<b>Base</b>			
<b>Height</b>	2.2" (5.6 cm)	2.2" (5.6 cm)	2.2" (5.6 cm)
<b>Length</b>	32" (81 cm)	32" (81 cm)	32" (81 cm)
<b>Width</b>	21.25" (54 cm)	21.25" (54 cm)	21.25" (54 cm)
<b>Caster diameter</b>	1.46" (3.7 cm)	1.46" (3.7 cm)	1.46" (3.7 cm)

### Storage compartment

Part number 3151-050-001

	Single-top without vanity	Single-top with vanity	Split-top with vanity
<b>Weight capacity</b>	40 lb. (18 kg)	40 lb. (18 kg)	40 lb. (18 kg)

### Total unit with storage compartment

	Single-top without vanity	Single-top with vanity	Split-top with vanity
<b>Weight</b>	68 lb. (31 kg)	75 lb. (43 kg)	86 lb. (39 kg)
<b>Length</b>	42.75" (109 cm)	42.75" (109 cm)	42.75" (109 cm)
<b>Width</b>	17.5" (44.5 cm)	17.5" (44.5 cm)	17.5" (44.5 cm)

### Warranty

One year parts and labor

Mkt lit-95-02022005 Rev G  
Copyright © 2018 Stryker  
Printed in U.S.A.

3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770

[stryker.com](http://stryker.com)

\*25 lb. safe working load  
Furniture is manufactured exclusively for Stryker Medical.  
Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Dura-Edge, Tru-Fit, Stryker. All other trademarks are trademarks of their respective owners or holder.