

Bedside Stands



Contemporary

Concealed casters enhance aesthetics and help protect casters from breaking

Removable drawers and plastic drawer liners are extremely durable and make the stands easy to clean

Steel drawer frames and hinges come with a lifetime warranty

Specialized thermoform process creates a seam-free finish on top and front panels

Vented back allows moisture and odors within to air out

Raised edges on top prevents spills and facilitates cleaning

Designs for a Healing Environment

Durable yet elegant, Stryker bedside stands are designed to meet the demanding needs of today's clinical environment. With five unique styles and a variety of laminate and hardware options, these bedside stands provide design latitude to coordinate your patient suites.



Legend™



Premier™



Shaker



Traditional

Bedside Stands

stryker[®]



Medical

3800 E. Centre Ave.
Portage, MI 49002
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Legend Standard Features

- Concealed 1.5" (4 cm) dual-wheel casters
- Finished edges
- Plastic drawers and drawer frames
- Polished chrome or brass pull handles
- Solid T-band edgemold
- Vented back

Premier Standard Features

- Concealed 1.5" (4 cm) dual-wheel casters
- Finished edges
- Removable drawers and plastic drawer liners
- Polished chrome or brass pull handles
- Raised spill containment top
- Vented back
- Steel drawer frames and glides

Contemporary Standard Features

- Concealed 1.5" (4 cm) dual-wheel casters
- Finished edges
- Removable drawers and plastic drawer liners
- Vented back
- Steel drawer frames and glides
- Thermoform drawer fronts and top
- Stainless steel pull handles

Shaker Standard Features

- Concealed 1.5" (4 cm) twin-wheel casters
- Finished edges
- Removable drawers and plastic drawer liners
- Vented back
- Steel drawer frames and glides
- Thermoform drawer fronts and top
- Round brushed stainless or antique brass knobs
- Recessed drawer and door panels

Traditional Standard Features

- Concealed 1.5" (4 cm) twin-wheel casters
- Finished edges
- Removable drawers and plastic drawer liners
- Vented back
- Steel drawer frames and glides
- Thermoform drawer fronts and top
- Polished chrome or brass pull handles
- Raised drawer and door panels

Specifications

• Style	Legend	Premier	Contemporary	Shaker	Traditional
	Part #4400-500/503-000	Part #4400-600/603-000	Part #4400-800/803-000	Part #4400-840/843-000	Part #4400-880/883-000
	1 Drawer with 1 Door/3 Drawer	1 Drawer with 1 Door/3 Drawer	1 Drawer with 1 Door/3 Drawer	1 Drawer with 1 Door/3 Drawer	1 Drawer with 1 Door/3 Drawer
• Overall Width	18.5" (47 cm)	19.75" (50 cm)	20" (51 cm)	19.25" (49 cm)	19.25" (49 cm)
• Overall Height	30.5" (77 cm)	30.5" (77 cm)	30.5" (77 cm)	30.5" (77 cm)	30.5" (77 cm)
• Overall Depth	18" (46 cm)	19.75" (50 cm)	20" (51 cm)	18.5" (47 cm)	18.5" (47 cm)
• Weight	80 lbs (36 kg)	80 lbs (36 kg)	80 lbs (36 kg)	80 lbs (36 kg)	80 lbs (36 kg)

- BIFMA modified static load testing up to 250 lbs.

Stryker reserves the right to change specifications without notice.
Furniture is manufactured for Stryker Medical.

Optional Features

- Three-drawer or drawer/door configuration
- Spill containment top (excluding Legend)
- Continuous pull handles (Premier and Contemporary models only)
- Key or slam lock on top drawer
- Choose from seven standard laminates

Warranty

- One year parts and labor
- Lifetime warranty on steel drawer frames and hinges

Custom Capabilities

- Custom capabilities available - please contact your Stryker sales representative

Copyright © 2005 Stryker
LSS-457-13E
PC/SCG-7.5M-0907
Mkt Lit-96-02022005 Rev E